


Improving Resolution for the KATRIN Detector Prototype

By Tabatha Spencer

Ursinus College, Collegeville PA
CENPA, University of Washington

With special thanks to Keith Rielage, Peter Doe,
Greg Harper, Sean McGee, Lesley Reece and
the CENPA EWI Group



Basic Outline


- KATRIN
 - Purpose
 - Schematic
- Equipment
 - Silicon Detector
 - Electron Gun
 - Vacuum System
- Resolution Techniques
 - Dead Layer Determination
 - Noise reduction
- Cooling Designs
 - Attempts
 - Future plans



KATRIN


(Karlsruhe Tritium Neutrino Experiment)

- Next generation tritium beta-decay neutrino experiment
- Direct measurement of neutrino mass
- Collaboration in Karlsruhe, Germany
- Projected start date in 2006



High Resolution


tritium β -decay and the neutrino rest mass



superallowed


half life : $t_{1/2} = 12.32 \text{ a}$

β end point energy : $E_0 = 18.57 \text{ keV}$




Silicon Detectors


- Testing resolution refining techniques on less sophisticated technology first
- Initial - Silicon Pin Diode
- Intermediate – Silicon Drift Diode
- Final - Segmented Pin Diode




Electron Gun and Vacuum System




Initial Measurements



8-Bit Bias





Scale Calibration



Dead Layer Measurement


- Small amount of energy lost on surface of the detector
- Correction needed
- Peak energy compared at multiple angles



Dead Layer Results

$$\Delta E = dE/dx * \Delta x_w * (1/\cos\alpha - 1)$$

- Can solve for the dead layer thickness
- $\Delta x_w = 123 \text{ nm}$
- Agrees well within error of previous measurement




Next Step: Cooling

- Try to significantly reduce amount of noise
- Clear precise peaks required
- Operating at low temperatures is a good way to do this
- Inexpensive and available devices:
 - Chiller
 - Peltier Device



Cooling Plan



Repeated Failure and Death



Current and Future Cooling Plans

- Scrapped Peltier plans till replacements available
- Go straight from cold rod to detector
- Look for better flexible connectors than braid
 - Multiple braids
 - Copper foil
 - Other possibilities?

THE END!!